

Chip Array Ball Grid Array

ZiLOG offers a 144-lead chip array ball grid array package (Figure 45).

- 1. Mark Permanency 3X soak into trichlorethane 1.1.1
1 minute duration each soak
Mech. brush after each soak
- 2. Hermeticity 5 X 10E-8 CC/SEC
MIL-STD-883C Method 1014.8 Condition B

► **Note:** The package dimensions are given in inches. To convert inches to millimeters, multiply the dimension by 25.4.




Figure 45. 144-Lead Chip Array Ball Grid Array